

For Research Use Only

# GZMH Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 68931-2-PBS

## Basic Information

Catalog Number:

68931-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11741

GenBank Accession Number:

BC027974

GeneID (NCBI):

2999

UNIPROT ID:

P20718

Full Name:

granzyme H (cathepsin G-like 2, protein h-CCPX)

Calculated MW:

246 aa, 27 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G9A4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50355-1, GZMH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68931-1-PBS. Detection antibody: 68931-2-PBS. Standard:Ag11741. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50355-2, GZMH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68931-2-PBS. Detection antibody: 68931-3-PBS. Standard:Ag11741. Range: 6.25-100 ng/mL.