

For Research Use Only

# CLEC7A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68934-1-PBS



## Basic Information

Catalog Number:

68934-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19019

GenBank Accession Number:

BC013385

GeneID (NCBI):

64581

UNIPROT ID:

Q9BXN2

Full Name:

C-type lectin domain family 7,  
member A

Calculated MW:

247 aa, 28 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A1A8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50361-1, CLEC7A/Dectin-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68934-1-PBS. Detection antibody: 68934-2-PBS. Standard:Ag19019. Range: 0.098-50 ng/mL.