

For Research Use Only

# CKAP4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68938-1-PBS



## Basic Information

Catalog Number: 68938-1-PBS	GenBank Accession Number: BC082972	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 10970	CloneNo.: 1E4E3
Source: Mouse	UNIPROT ID: Q07065	
Isotype: IgG1	Full Name: cytoskeleton-associated protein 4	
Immunogen Catalog Number: AG10022	Calculated MW: 602 aa, 66 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50368-1, CKAP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68938-1-PBS. Detection antibody: 68938-2-PBS. Standard:Ag10022. Range: 1.563-100 ng/mL.