

For Research Use Only

Cytokeratin 5/6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68944-2-PBS



Basic Information

Catalog Number: 68944-2-PBS	GenBank Accession Number: BC034535	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 3854	CloneNo.: 2E10E4
Source: Mouse	UNIPROT ID: P04259	
Isotype: IgG2b	Full Name: keratin 6B	
Immunogen Catalog Number: AG11241	Calculated MW: 564 aa, 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50379-1, KRT6B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68944-1-PBS. Detection antibody: 68944-2-PBS. Standard: Ag11241. Range: 0.781-50 ng/mL