

For Research Use Only

PPIF-Specific Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68956-5-PBS



Basic Information

Catalog Number:

68956-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG13151

GenBank Accession Number:

BC005020

GeneID (NCBI):

10105

UNIPROT ID:

P30405

Full Name:

peptidylprolyl isomerase F

Calculated MW:

22 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B7B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50406-3, PPIF-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68956-5-PBS. Detection antibody: 68956-3-PBS. Standard: Ag13151. Range: 1.563-100 ng/mL.