

For Research Use Only

# TNFRSF14 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68961-1-PBS



## Basic Information

Catalog Number: 68961-1-PBS	GenBank Accession Number: BC002794	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 8764	CloneNo.: 2A4A9
Source: Mouse	UNIPROT ID: Q92956	
Isotype: IgG1	Full Name: tumor necrosis factor receptor superfamily, member 14 (herpesvirus entry mediator)	
Immunogen Catalog Number: AG28731	Calculated MW: 30 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50404-5, TNFRSF14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68961-1-PBS. Detection antibody: 68961-2-PBS. Standard:Ag28731. Range: 0.098-3.125 ng/mL