

For Research Use Only

# CD71 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68962-1-PBS



## Basic Information

Catalog Number:

68962-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

7037

ENSEMBL Gene ID:

ENSG00000072274

UNIPROT ID:

P02786

Full Name:

transferrin receptor (p90, CD71)

Calculated MW:

85 kDa

Purification Method:

Protein G purification

CloneNo.:

1B10D8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50414-1, CD71 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68962-1-PBS. Detection antibody: 66180-2-PBS. Standard:Eg0339. Range: 1.563-50 ng/mL



Cytometric bead array standard curve of MP50414-2, CD71 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68962-1-PBS. Detection antibody: 68962-2-PBS. Standard:Eg0339. Range: 0.781-100 ng/mL