

For Research Use Only

GPR64 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68963-2-PBS



Basic Information

Catalog Number:

68963-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31274

GenBank Accession Number:

BC099901

GeneID (NCBI):

10149

UNIPROT ID:

Q8IZP9

Full Name:

G protein-coupled receptor 64

Calculated MW:

1017 aa, 112 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D10C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50415-1, GPR64 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68963-1-PBS. Detection antibody: 68963-2-PBS. Standard:Ag31274. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50415-2, GPR64 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68963-3-PBS. Detection antibody: 68963-2-PBS. Standard:Ag31274. Range: 0.195-100 ng/mL