

For Research Use Only

# ELOVL3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68964-2-PBS



## Basic Information

Catalog Number:

68964-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24318

GenBank Accession Number:

BC034344

GeneID (NCBI):

83401

UNIPROT ID:

Q9HB03

Full Name:

elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 3

Calculated MW:

32 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D4G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50418-1, ELOVL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68964-1-PBS. Detection antibody: 68964-2-PBS. Standard:Ag24318. Range: 0.098-100 ng/mL.