

For Research Use Only

ELOVL3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68964-5-PBS



Basic Information

Catalog Number:

68964-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24318

GenBank Accession Number:

BC034344

GeneID (NCBI):

83401

UNIPROT ID:

Q9HB03

Full Name:

elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 3

Calculated MW:

32 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E12A12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50418-3, ELOVL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68964-5-PBS. Detection antibody: 68964-6-PBS. Standard:Ag24318. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50418-4, ELOVL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68964-5-PBS. Detection antibody: 68964-7-PBS. Standard:Ag24318. Range: 0.781-100 ng/mL.