

For Research Use Only

Integrin alpha-6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68968-3-PBS



Basic Information

Catalog Number:

68968-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG16335

GenBank Accession Number:

BC136456

GeneID (NCBI):

3655

UNIPROT ID:

P23229

Full Name:

integrin, alpha 6

Calculated MW:

1073 aa, 119 kDa

Purification Method:

Protein A purification

CloneNo.:

2B11E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50396-5, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68968-2-PBS. Detection antibody: 68968-3-PBS. Standard:Ag16335. Range: 0.781-100 ng/mL.