

For Research Use Only

ADAMTS17 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68973-3-PBS



Basic Information

Catalog Number:

68973-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG32013

GenBank Accession Number:

NM_139057.3

GeneID (NCBI):

170691

UNIPROT ID:

Q8TE56

Full Name:

ADAM metalloproteinase with
thrombospondin type 1 motif, 17

Calculated MW:

121 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A4D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50434-2, ADAMTS17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68973-1-PBS. Detection antibody: 68973-3-PBS. Standard:Ag32013. Range: 0.781-100 ng/mL.