

For Research Use Only

MLXIP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68977-2-PBS



Basic Information

Catalog Number: 68977-2-PBS	GenBank Accession Number: BC039704	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 22877	CloneNo.: 2G12F2
Isotype: IgG1	UNIPROT ID: Q9HAP2	
Immunogen Catalog Number: AG4519	Full Name: MLX interacting protein	
	Calculated MW: 919 aa, 100 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51488-1, MLXIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68977-2-PBS. Detection antibody: 68977-3-PBS.



Cytometric bead array standard curve of MP51488-1, MLXIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68977-2-PBS. Detection antibody: 68977-3-PBS. Standard: Ag4519. Range: 0.781-100 ng/mL.