

For Research Use Only

Osteopontin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68983-5-PBS



Basic Information

Catalog Number:

68983-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19216

GenBank Accession Number:

BC007016

GeneID (NCBI):

6696

UNIPROT ID:

P10451-5

Full Name:

secreted phosphoprotein 1

Calculated MW:

314 aa, 35 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H7E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array of Human urine using MP50455-3, Osteopontin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68983-4-PBS. Detection antibody: 68983-5-PBS. Urine of four individual healthy human donors was analyzed using MP50455-3.



Cytometric bead array standard curve of MP50455-3, OPN, SPP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68983-4-PBS. Detection antibody: 68983-5-PBS. Standard: Ag19216. Range: 3.125-100 ng/mL.