

For Research Use Only

CEP170 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68987-3-PBS



Basic Information

Catalog Number:

68987-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG26318

GenBank Accession Number:

BC140794

GeneID (NCBI):

9859

UNIPROT ID:

Q55W79

Full Name:

centrosomal protein 170kDa

Calculated MW:

1486 aa, 165 kDa

Purification Method:

Euglobulin precipitation

CloneNo.:

1D7A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50467-2, CEP170 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68987-1-PBS. Detection antibody: 68987-3-PBS. Standard:Ag26318. Range: 0.781-100 ng/mL.