

For Research Use Only

BIM, BCL2L11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68988-1-PBS



Basic Information

Catalog Number:

68988-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG20546

GenBank Accession Number:

BC033694

GeneID (NCBI):

10018

UNIPROT ID:

O43521

Full Name:

BCL2-like 11 (apoptosis facilitator)

Calculated MW:

198 aa, 22 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

10A10E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50468-1, BIM, BCL2L11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68988-1-PBS. Detection antibody: 68988-2-PBS. Standard:Ag20546. Range: 0.098-12.5 ng/mL