

For Research Use Only

MRPL11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68992-1-PBS



Basic Information

Catalog Number:

68992-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8114

GenBank Accession Number:

BC005002

GeneID (NCBI):

65003

UNIPROT ID:

Q9Y3B7

Full Name:

mitochondrial ribosomal protein L11

Calculated MW:

21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F4A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50482-1, MRPL11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68992-1-PBS. Detection antibody: 68992-2-PBS. Standard:Ag8114. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50482-3, MRPL11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68992-1-PBS. Detection antibody: 68992-3-PBS. Standard:Ag8114. Range: 0.781-100 ng/mL