

For Research Use Only

Jagged1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68997-1-PBS



Basic Information

Catalog Number:

68997-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28216

GenBank Accession Number:

BC098393

GeneID (NCBI):

182

UNIPROT ID:

P78504

Full Name:

jagged 1 (Alagille syndrome)

Calculated MW:

1218 aa, 134 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F8E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50489-1, Jagged1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68997-1-PBS. Detection antibody: 68997-2-PBS. Standard:Ag28216. Range: 0.195-100 ng/mL.