

For Research Use Only

IFN alpha 2B Monoclonal antibody

Catalog Number: 69016-1-Ig



Basic Information

Catalog Number:

69016-1-Ig

Size:

N/A

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1072

GenBank Accession Number:

GeneID (NCBI):

3440

Full Name:

interferon, alpha 2

Purification Method:

Protein G purification

CloneNo.:

1E5C4

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Lyophilized antibodies are stable for 1 year from the date of receipt if stored between (-20°C) and (-80°C). Upon reconstitution we recommend that the solution can be stored at(4°C) for short term or at(-20°C) to (-80°C) for long term. Repeated freeze thaw cycles should be avoided with reconstituted products.

Storage Buffer:

Sterile PBS.

Aliquoting is unnecessary for -20° C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

