

For Research Use Only

IL-10 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 69018-2-PBS



Basic Information

Catalog Number:

69018-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

HZ-1145

GenBank Accession Number:

GeneID (NCBI):

3586

Full Name:

interleukin 10

Purification Method:

Protein A purification

CloneNo.:

7B8G3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50035-1, IL-10
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 69018-1-PBS. Detection
antibody: 69018-2-PBS. Standard: HZ-1145. Range:
15.6-1000 pg/mL.