

For Research Use Only

# rat anti HA Monoclonal antibody (7C9)



[www.ptgcn.com](http://www.ptgcn.com)

Catalog Number: 7c9

28 Publications

## Basic Information

**Catalog Number:**  
7c9

**Clone No.:**  
7C9

**Applications:**  
WB, IF, IP

**Host/IsoType:**  
Rat/IgG2a

**Conjugate:**  
Unconjugated

**Type:**  
Primary Antibody

**Class:**  
Monoclonal

**RRID:**  
AB\_2631399

**Purification Method:**  
Protein G purified

## Description

Rat monoclonal [7C9] to HA-tag

## Specificity/Target

HA-tag (YPYDVPDYA) at the N-terminus, C-terminus, or internal site of fusion protein.

## Recommended Dilution

Western Blot: 1:1,000  
Immunofluorescence: 1:1,000  
Immunoprecipitation: 5 µg/ IP

## Storage

**Storage:**  
Shipped on ice. Store at -20°C/-4°F. Stable for 1 year after shipment.

**Storage Buffer:**  
PBS, preservative: 0.09 % Sodium azide


For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)  
W: [www.ptgcn.com](http://www.ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Western blot analysis of cell extract from HEK293T cells expressing GFP-HA (29.9 kDa)



Immunofluorescence of HeLa cells transiently expressing HA-TOMM70