

For Research Use Only

# mouse IgG Kappa chain Recombinant antibody



Catalog Number: 80015-1-RR

## Basic Information

Catalog Number:

80015-1-RR

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

155 kDa

Purification Method:

Protein A purification

CloneNo.:

3D9

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HEK-293 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HEK-293 cells were subjected to SDS-PAGE followed by western blot with eight different monoclonal antibodies and subsequently detected by recombinant rabbit anti-mouse kappa chain antibody 80015-1-Ig at 1:10000 and then probed with HRP-goat anti-rabbit secondary antibody.