

For Research Use Only

Phospho-p70 S6 kinase (Thr389) Recombinant monoclonal antibody

Catalog Number: 80021-1-RR

1 Publications



Basic Information

Catalog Number:

80021-1-RR

Concentration:

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003161

GeneID (NCBI):

6198

UNIPROT ID:

P23443

Full Name:

ribosomal protein S6 kinase, 70kDa,
polypeptide 1

Purification Method:

Protein A purification

CloneNo.:

4C8

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xue Jiao	40367705	Vet Microbiol	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

