

For Research Use Only

# STAT2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 80061-2-PBS



## Basic Information

Catalog Number:

80061-2-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10168

GenBank Accession Number:

BC051284

GeneID (NCBI):

6773

UNIPROT ID:

P52630

Full Name:

signal transducer and activator of transcription 2, 113kDa

Calculated MW:

851 aa, 98 kDa

Purification Method:

Protein A purification

CloneNo.:

242685B9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01853-2, STAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80061-2-PBS. Detection antibody: 80061-4-PBS. Standard: Ag10168. Range: 1.563-200 ng/mL.