

For Research Use Only

Phospho-TAU (Thr181) Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 80073-5-PBS



Basic Information

Catalog Number:

80073-5-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC000558

GeneID (NCBI):

4137

UNIPROT ID:

P10636

Full Name:

microtubule-associated protein tau

Calculated MW:

37-46, 79-81 kDa

Purification Method:

Protein A purification

CloneNo.:

230429H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00210-1, TAU Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80073-5-PBS. Detection antibody: 80073-2-PBS. Standard: SY01799P. Range: 0.156-2.5 ng/mL



Cytometric bead array standard curve of MP00210-2, TAU Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80073-4-PBS. Detection antibody: 80073-5-PBS. Standard: SY01799P. Range: 0.156-2.5 ng/mL



Cytometric bead array standard curve of MP00210-3, TAU Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80073-3-PBS. Detection antibody: 80073-5-PBS. Standard: SY01799P. Range: 0.156-2.5 ng/mL