

For Research Use Only

Phospho-AMPK Alpha (Thr172) Recombinant monoclonal antibody

Catalog Number: 80209-6-RR

2 Publications



Basic Information

Catalog Number:

80209-6-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC048980

GeneID (NCBI):

5562

UNIPROT ID:

Q13131

Full Name:

protein kinase, AMP-activated, alpha
1 catalytic subunit

Observed MW:

64 kDa

Purification Method:

Protein A purification

CloneNo.:

242923D4

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB: H₂O₂ treated C2C12 cells, λ phosphatase treated
HEK-293 cells

Background Information

AMPK is a serine/threonine protein kinase complex consisting of a catalytic α -subunit (α 1 and α 2), a scaffolding β -subunit (β 1 and β 2), and a regulatory γ -subunit (γ 1, γ 2, and γ 3). Ubiquitous expression of AMPK α 1-, β 1-, and γ 1-subunits in many tissues makes the α 1 β 1 γ 1 complex a reference for AMPK assays to identify AMPK activators. AMPK is generally quiescent under normal conditions but is activated in response to signals and stresses that increase the AMP/ATP ratio, such as hypoglycemia, strenuous exercise, anoxia, and ischemia. An increase in the ratio of AMP/ATP activates AMPK by several mechanisms, including direct allosteric activation and covalent modification in which an AMP-dependent AMPK kinase (AMPKK) phosphorylates the α -subunit on Thr172. Once activated, AMPK switches on catabolic pathways that generate ATP, while switching off ATP-consuming processes (e.g., biosynthesis, cell growth, and proliferation), and in doing so acts as an "energy gauge". (PMID: 27034026, PMID: 21980456, PMID: 27600021)

This antibody can recognize the phosphorylation sites of Thr183 in AMPK Alpha 1 and Thr172 in AMPK Alpha 2.

Notable Publications

Author	Pubmed ID	Journal	Application
Cheng Yu	40500724	Mol Med	WB
Yongmei Zhang	40450199	BMC Cardiovasc Disord	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Non-treated C2C12 cells and H2O2 treated C2C12 cells were subjected to SDS PAGE followed by western blot with 80209-6-RR (Phospho-AMPK Alpha (Thr172) antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with GAPDH (60004-1-Ig) antibody as a loading control.



Non-treated, phosphatase inhibitor treated and λ phosphatase treated HEK-293 cells were subjected to SDS PAGE followed by western blot with 80209-6-RR (Phospho-AMPK Alpha (Thr172) antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Beta Actin (66009-1-Ig) antibody as a loading control.