

For Research Use Only

TTYH1 Recombinant antibody

Catalog Number: 80263-1-RR



Basic Information

Catalog Number:

80263-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25290

GenBank Accession Number:

BC011347

GeneID (NCBI):

57348

UNIPROT ID:

Q9H313

Full Name:

tweety homolog 1 (Drosophila)

Calculated MW:

49-50 kDa

Purification Method:

Protein A purification

CloneNo.:

8G21

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, rabbit, pig

Positive Controls:

WB : pig cerebellum tissue, human testis tissue, rabbit cerebellum tissue, rabbit brain tissue

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80263-1-RR (TTYH1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.