

For Research Use Only

# XRCC5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 80325-4-PBS

## Basic Information

Catalog Number: 80325-4-PBS	GenBank Accession Number: BC019027	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7520	CloneNo.: 240301B10
Source: Rabbit	UNIPROT ID: P13010	
Isotype: IgG	Full Name: X-ray repair complementing defective repair in Chinese hamster cells 5 (double-strand-break rejoining)	
Immunogen Catalog Number: AG9454	Calculated MW: 732 aa, 83 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00384-1, XRCC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80325-5-PBS. Detection antibody: 80325-4-PBS. Standard: Ag9454. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00384-3, XRCC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80325-3-PBS. Detection antibody: 80325-4-PBS. Standard: Ag9454. Range: 0.625-80 ng/mL