

For Research Use Only

XRCC5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 80325-5-PBS

Basic Information

Catalog Number:

80325-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9454

GenBank Accession Number:

BC019027

GeneID (NCBI):

7520

UNIPROT ID:

P13010

Full Name:

X-ray repair complementing defective repair in Chinese hamster cells 5 (double-strand-break rejoining)

Calculated MW:

732 aa, 83 kDa

Purification Method:

Protein A purification

CloneNo.:

240301F8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00384-1, XRCC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80325-5-PBS. Detection antibody: 80325-4-PBS. Standard: Ag9454. Range: 0.625-80 ng/mL