

For Research Use Only

Perilipin-2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 80362-4-PBS



Basic Information

Catalog Number:

80362-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7539

GenBank Accession Number:

BC005127

GeneID (NCBI):

123

UNIPROT ID:

Q99541

Full Name:

adipose differentiation-related protein

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230274E5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00071-1, Perilipin-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80362-3-PBS. Detection antibody: 80362-4-PBS. Standard: Ag7539. Range: 0.156-10 ng/mL