

For Research Use Only

NLRP3 Recombinant antibody

Catalog Number: 80422-2-RR



Basic Information

Catalog Number:	80422-2-RR	GenBank Accession Number:	NM_001079821	Purification Method:	Protein A purification
Size:		GeneID (NCBI):	114548	CloneNo.:	230052F12
Source:	Rabbit	Full Name:	NLR family, pyrin domain containing 3		
Isotype:	IgG	Calculated MW:	118 kDa		
Immunogen Catalog Number:	AG26289				

Applications

Tested Applications:
ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

