

For Research Use Only

E-cadherin Recombinant antibody, PBS Only (Detector)



Catalog Number: 80541-2-PBS

Basic Information

Catalog Number:

80541-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14973

GenBank Accession Number:

BC141838

GeneID (NCBI):

999

UNIPROT ID:

P12830

Full Name:

cadherin 1, type 1, E-cadherin
(epithelial)

Calculated MW:

882 aa, 97 kDa

Purification Method:

Protein A purification

CloneNo.:

240159A5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00187-1, E-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80541-5-PBS. Detection antibody: 80541-2-PBS. Standard: Ag14973. Range: 0.625-20 ng/mL.