

For Research Use Only

E-cadherin Recombinant antibody, PBS Only (Detector)



Catalog Number: 80541-6-PBS

Basic Information

Catalog Number: 80541-6-PBS	GenBank Accession Number: BC141838	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 999	CloneNo.: 240159E2
Source: Rabbit	UNIPROT ID: P12830	
Isotype: IgG	Full Name: cadherin 1, type 1, E-cadherin (epithelial)	
Immunogen Catalog Number: AG14973	Calculated MW: 882 aa, 97 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00187-2, E-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80541-4-PBS. Detection antibody: 80541-6-PBS. Standard: Ag14973. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00187-3, E-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80541-3-PBS. Detection antibody: 80541-6-PBS. Standard: Ag14973. Range: 1.25-20 ng/mL.