

For Research Use Only

# Tubulin-beta Recombinant antibody, PBS Only (Capture)



Catalog Number: 80680-1-PBS

## Basic Information

Catalog Number: 80680-1-PBS	GenBank Accession Number: BC000748	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10381	CloneNo.: 2E8
Source: Rabbit	UNIPROT ID: Q13509	
Isotype: IgG	Full Name: tubulin, beta 3	
Immunogen Catalog Number: AG0117	Calculated MW: 450 aa, 50 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sample test of MP00550-1, Tubulin-beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80680-1-PBS. Detection antibody: 80676-1-PBS. Cell lysate(30 µ g/well): HEK293, Jurkat, K-562, HSC-T6, NIH/3T3, RAW264.7.



Cytometric bead array standard curve of MP00550-1, Tubulin-beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80680-1-PBS. Detection antibody: 80676-1-PBS. Standard:Ag0117. Range: 0.098-100 ng/mL.