

For Research Use Only

ANGPTL4 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 80738-3-PBS



Basic Information

Catalog Number:

80738-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5329

GenBank Accession Number:

NM_139314.2

GeneID (NCBI):

51129

UNIPROT ID:

Q9BY76-1

Full Name:

angiopoietin-like 4

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

252229E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ANGPTL4 concentration was determined to be 17.0 ng/mL in human serum.



Cytometric bead array standard curve of MP02878-1, ANGPTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80738-3-PBS. Detection antibody: 80738-2-PBS. Standard: Eg5329. Range: 0.469-60 ng/mL.