

For Research Use Only



Tubulin-Alpha Recombinant antibody, PBS Only (Capture)

Catalog Number: 80740-1-PBS

Basic Information

Catalog Number: 80740-1-PBS	GenBank Accession Number: BC009314	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10376	CloneNo.: 3H19
Source: Rabbit	UNIPROT ID: P68363	
Isotype: IgG	Full Name: tubulin, alpha 1b	
Immunogen Catalog Number: AG18034	Calculated MW: 50 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00551-1, Tubulin-Alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80740-1-PBS. Detection antibody: 80762-1-PBS. Standard:Ag18034. Range: 0.098-50 ng/mL.