

For Research Use Only



# DYKDDDDK tag Recombinant antibody, PBS Only (Binds to FLAG® tag epitope)

Catalog Number: **80745-1-PBS**

## Basic Information

Catalog Number:

80745-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2329

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

5E11

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

recombinant protein

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



TDP-43-flag Transfected HEK-293T cells were subjected to SDS PAGE followed by western blot with 80745-1-RR (DYKDDDDK tag antibody) at dilution of 1:50000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80745-1-PBS in a different storage buffer formulation.