

For Research Use Only



DYKDDDDK tag Recombinant antibody (Binds to FLAG® tag epitope)

Catalog Number: **80745-1-RR**

Basic Information

Catalog Number:

80745-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2329

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

5E11

Recommended Dilutions:

WB 1:20000-1:100000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

recombinant protein

Positive Controls:

WB : TDP-43-flag Transfected HEK-293T cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



TDP-43-flag Transfected HEK-293T cells were subjected to SDS PAGE followed by western blot with 80745-1-RR (DYKDDDDK tag antibody) at dilution of 1:50000 incubated at room temperature for 1.5 hours.