

For Research Use Only

# beta IG-H3/TGFBI Recombinant antibody, PBS Only (Capture)

Catalog Number: 80805-3-PBS



## Basic Information

Catalog Number:

80805-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1233

GenBank Accession Number:

NM\_000358.3

GeneID (NCBI):

7045

ENSEMBL Gene ID:

ENSG00000120708

UNIPROT ID:

Q15582

Full Name:

transforming growth factor, beta-induced, 68kDa

Calculated MW:

75kDa

Purification Method:

Protein A purification

CloneNo.:

243108H4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean TGFBI/BIGH3 concentration was determined to be 3586.8 ng/mL in human serum.



Cytometric bead array standard curve of MP02279-1, TGFBI/BIGH3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80805-3-PBS. Detection antibody: 80805-2-PBS. Standard: Eg1233. Range: 0.781-100 ng/mL.