

For Research Use Only

GSDMD Recombinant antibody, PBS Only (Detector)

Catalog Number: 80835-6-PBS



Basic Information

Catalog Number:

80835-6-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC008904

GeneID (NCBI):

79792

UNIPROT ID:

P57764

Full Name:

gasdermin D

Calculated MW:

484 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

242459H7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02194-1, Human GSDMD Recombinant Matched Antibody Pair - PBS only. 80835-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag35247. 80835-6-PBS was HRP conjugated as the detection antibody. Range: 1.56-50 ng/mL



The mean GSDMD concentration was determined to be 250.44 ng/mL in THP-1 cell extract based on a 1.50 mg/mL extract load.