

For Research Use Only

LAG3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 80855-4-PBS



Basic Information

Catalog Number:

80855-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002286

GeneID (NCBI):

3902

UNIPROT ID:

P18627

Full Name:

lymphocyte-activation gene 3

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

243086C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02008-1, LAG-3/CD223 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80855-4-PBS. Detection antibody: 80855-3-PBS. Standard: Eg31550. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP02008-2, LAG-3/CD223 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80855-4-PBS. Detection antibody: 80855-2-PBS. Standard: Eg31550. Range: 0.781-100 ng/mL