

For Research Use Only

Fibronectin Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 80911-10-PBS



Basic Information

Catalog Number:

80911-10-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5117

GenBank Accession Number:

NM_212482.4

GeneID (NCBI):

2335

UNIPROT ID:

P02751-15

Full Name:

fibronectin 1

Calculated MW:

272 kDa

Purification Method:

Protein A purification

CloneNo.:

251989E9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




HepG2 (human hepatocellular carcinoma cells) were cultured in DMEM supplemented with 10% fetal bovine serum, 2.5 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. An aliquot of the cell culture supernate was removed, assayed for human FN1, and measured 2,102.36 ng/mL.




Saliva of six individual healthy human donors was measured. The human FN1 concentration of detected samples was determined to be 38.99 ng/mL with a range of 17.55 - 64.99 ng/mL.




Urine of five individual healthy human donors was measured. The human FN1 concentration of detected samples was determined to be 18.35 ng/mL with a range of 8.2 - 25.2 ng/mL.




Human milk of six individual healthy human donors was measured. The human FN1 concentration of detected samples was determined to be 224.52 ng/mL with a range of 55.37 - 371.75 ng/mL.



The mean FN1 concentration was determined to be 136.80 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Plasma of six individual healthy human donors was measured. The human FN1 concentration of detected samples was determined to be 70,005.24 ng/mL with a range of 37,468.04 - 107,177.08 ng/mL.



Sandwich ELISA standard curve of MP02923-1, Human FN1 Recombinant Matched Antibody Pair - PBS only. 80911-10-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5117. 80911-9-PBS was HRP conjugated as the detection antibody. Range: 78.1-2500 pg/mL.