

For Research Use Only

Fibronectin Recombinant monoclonal antibody, PBS Only

Catalog Number: 80911-4-PBS



Basic Information

Catalog Number:

80911-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5117

GenBank Accession Number:

NM_212482.4

GeneID (NCBI):

2335

UNIPROT ID:

P02751-15

Full Name:

fibronectin 1

Calculated MW:

272 kDa

Observed MW:

220-275 kDa

Purification Method:

Protein A purification

CloneNo.:

251989E12

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Fibronectin (fibronectin) is a kind of macromolecular glycoprotein widely existing in extracellular matrix (ECM), plasma and various cell surfaces, which plays a central role in key biological processes such as cell adhesion, migration, proliferation, differentiation, tissue repair and embryo development. Fibronectin is usually a dimer structure, which consists of two subunits connected by disulfide bonds. There are two main forms of Fibronectin: plasma Fibronectin(pFN) is synthesized by liver cells and exists in the blood, which is involved in the processes of coagulation and wound repair. Cellular Fibronectin(cFN) is synthesized by fibroblasts and other cells, deposited in extracellular matrix, and participates in the interaction between cells and matrix. Its molecular weight ranges from 220 to 275 kDa.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 80911-4-RR (Fibronectin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80911-4-PBS in a different storage buffer formulation.



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 80911-4-RR (Fibronectin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80911-4-PBS in a different storage buffer formulation.



rat plasma were subjected to SDS PAGE followed by western blot with 80911-4-RR (Fibronectin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80911-4-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed NIH/3T3 cells using Fibronectin antibody (80911-4-RR, Clone: 251989E12) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 80911-4-PBS in a different storage buffer formulation.