

For Research Use Only

FGFR4 Recombinant antibody, PBS Only



Catalog Number: 81069-1-PBS

Basic Information

Catalog Number: 81069-1-PBS	GenBank Accession Number: BC011847	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2264	CloneNo.: 5P12
Source: Rabbit	UNIPROT ID: P22455	
Isotype: IgG	Full Name: fibroblast growth factor receptor 4	
Immunogen Catalog Number: AG31013	Calculated MW: 88 kDa	
	Observed MW: 100-110 kDa	

Applications

Tested Applications:
Indirect ELISA, WB

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81069-1-RR (FGFR4 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81069-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 81069-1-RR against Human FGFR4 were performed. The affinity constant is 25 nM.