

For Research Use Only

FGFR4 Recombinant antibody

Catalog Number: 81069-1-RR



Basic Information

Catalog Number: 81069-1-RR	GenBank Accession Number: BC011847	Purification Method: Protein A purification
Size: 800 µg/ml	GeneID (NCBI): 2264	CloneNo.: 5P12
Source: Rabbit	UNIPROT ID: P22455	Recommended Dilutions: WB 1:2000-1:16000
Isotype: IgG	Full Name: fibroblast growth factor receptor 4	
Immunogen Catalog Number: AG31013	Calculated MW: 88 kDa	
	Observed MW: 100-110 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : MOLT-4 cells, Raji cells, Jurkat cells, HeLa cells, HepG2 cells, HT-29 cells, A549 cells
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81069-1-RR (FGFR4 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.