

For Research Use Only

FGFR4 Recombinant antibody

Catalog Number: 81070-1-RR



Basic Information

Catalog Number:

81070-1-RR

Size:

2292 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31013

GenBank Accession Number:

BC011847

GeneID (NCBI):

2264

UNIPROT ID:

P22455

Full Name:

fibroblast growth factor receptor 4

Calculated MW:

88 kDa

Observed MW:

100-110 kDa

Purification Method:

Protein A purification

CloneNo.:

2M22

Recommended Dilutions:

WB 1:4000-1:40000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : MOLT-4 cells, Raji cells, Jurkat cells, HeLa cells, HepG2 cells, HT-29 cells, A549 cells

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81070-1-RR (FGFR4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.