

For Research Use Only

FASN Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:81079-4-PBS



Basic Information

Catalog Number: 81079-4-PBS	GenBank Accession Number: BC007909	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2194	CloneNo.: 260308D1
Isotype: IgG	UNIPROT ID: P49327	
Immunogen Catalog Number: AG0975	Full Name: fatty acid synthase	
	Calculated MW: 272 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FASN concentration was determined to be 17.3 ng/mL in HEK-293T cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03330-1, FASN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81079-4-PBS. Detection antibody: 81079-3-PBS. Standard: Ag0975. Range: 0.469-60 ng/mL.