

For Research Use Only

YAP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 81087-1-PBS



Basic Information

Catalog Number:

81087-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4510

GenBank Accession Number:

BC038235

GeneID (NCBI):

10413

UNIPROT ID:

P46937

Full Name:

Yes-associated protein 1, 65kDa

Calculated MW:

504 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

5H2

Affinity:

$K_D = 1.01 \times 10^{-10}$

$K_{off} = 4.24 \times 10^{-5}$

$K_{on} = 4.18 \times 10^5$

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00309-2, YAP1
Recombinant Matched Antibody Pair - PBS only.
Capture antibody: 81096-1-PBS. Detection
antibody: 81087-1-PBS. Standard: Ag4510. Range:
78.1-5000 pg/mL.