

For Research Use Only

YAP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 81096-1-PBS



Basic Information

Catalog Number: 81096-1-PBS	GenBank Accession Number: BC038235	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10413	CloneNo.: 1M10
Source: Rabbit	UNIPROT ID: P46937	Affinity: $K_D = 1.37 \times 10^{-10}$
Isotype: IgG	Full Name: Yes-associated protein 1, 65kDa	$K_{Off} = 1.08 \times 10^{-4}$
Immunogen Catalog Number: AG4510	Calculated MW: 504 aa, 54 kDa	$K_{On} = 7.87 \times 10^5$

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00309-1, YAP1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 81096-1-PBS. Detection antibody: 81099-1-PBS. Standard: Ag4510. Range: 78.1-5000 pg/mL



Sandwich ELISA standard curve of MP00309-2, YAP1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 81096-1-PBS. Detection antibody: 81087-1-PBS. Standard: Ag4510. Range: 78.1-5000 pg/mL