

For Research Use Only

AQP1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 81137-5-PBS



Basic Information

Catalog Number:

81137-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14093

GenBank Accession Number:

BC022486

GeneID (NCBI):

358

UNIPROT ID:

P29972

Full Name:

aquaporin 1 (Colton blood group)

Calculated MW:

269 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

251798C2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean AQP1 concentration was determined to be 1.74 ng/mL in mouse skeletal tissue extract based on a 4.1 mg/mL extract load and 89.24 ng/mL in mouse kidney tissue extract based on a 6.4 mg/mL extract load.



Plasma of six individual healthy human donors was measured. The human AQP1 concentration of detected samples was determined to be 435.6 pg/mL with a range of 126.0 - 853.9 pg/mL.



Sandwich ELISA standard curve of MP02804-1, Human AQP1 Recombinant Matched Antibody Pair - PBS only. 81137-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag14093. 81137-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL